

Apelin-12 (Human, Rat, Mouse, Bovine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-057-23	1 kit	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Apelin - Research Abstracts, Antibody Stainings, Kits Standard Curve, Sequences, etc.](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Assay Protocol \(PDF\)](#) - **No Sample Extraction Required***

*Serum and plasma do not require extraction. All other sample types may need extraction.

Sample Preparation [Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-057-23

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.07ng/ml

Range 0-100 ng/ml

Linear Range 0.07-0.79 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Apelin-12 (Human)	100
	Apelin-13 (Human)	100
	Apelin-13 (Rat)	100
	Apelin-36 (Human)	100
	Apelin-36 (Rat)	100
	Alpha-ANP (Human)	0
	ANP (25-56) (N-ANP)	0
	ADM-52 (Human)	0
	BNP-32 (Human)	0
	CNP-22	0
	Ghrelin (Human)	0
	ET-1 (Human)	0
	Bradykinin	0

