

Somatostatin (Human, Rat, Mouse, Porcine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-060-03	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Neuronostatin

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-060-03

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.13ng/ml

Range 0-25ng/ml

Linear Range 0.13-1.16 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Somatostatin	100
	Somatostatin-28	100
	Somatostatin-25	100
	[Des-Ala1]-Somatostatin	100
	Somatostatin (1-32)	0
	Pro-Somatostatin-32	0
	Somatostatin Analog, RC-160	0
	Somatostatin Analog (CTOP)	0
	Octreotide	0
	Insulin (Human)	0
	Glucagon (Human)	0
	Amylin (Human)	0
	Substance P	0
	VIP (Human, Porcine, Rat)	0
	NPY (Porcine)	0

Sequence Ala - Gly - Cys - Lys - Asn - Phe - Phe - Trp - Lys - Thr - Phe - Thr - Ser - Cys

Disulfide Bridge: Cys3 - Cys14

Standard
Curve Image

