

Neuronostatin-13 (Human, Porcine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-060-50	1 kit	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

Neuronostatin

Protocol	Assay Protocol (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	EK-060-50
Standard Size	96 wells
Family	Neuronostatin
Range	0-100 ng/ml
Linear Range	0.12-1.47 ng/ml

Thursday 20 December, 2012

11983462 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.