

# Vasopressin (AVP) [Arg8] (Human, Rat, Mouse, Ovine, Canine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-065-07	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

### Normal Plasma Peptide/Protein Levels

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-065-07

**Standard Size** 96 wells

**Sensitivity** Minimum Detectable Concentration = 0.04 ng/ml

**Range** 0-100 ng/ml


**Linear Range** 0.04-0.51 ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	[Arg8]-Vasopressin (Human, Rat, Mouse, Bovine, Ovine, Canine)	100
	[Arg8]-Vasotocin	78
	[Lys8]-Vasopressin (Human, Rat, Mouse, Ovine)	0
	Deamino-[D-Arg8]-Vasopressin	0
	Oxytocin	0
	LH-RH	0
	ACTH (Human)	0
	$\alpha$ -ANP(1-28) (Human, Canine)	0
	Met-Enkephalin	0
	PACAP38 (Human, Ovine, Rat)	0
	PACAP27-NH2 (Human, Ovine, Rat)	0

**Sequence** Cys - Tyr - Phe - Gln - Asn - Cys - Pro - Arg - Gly - NH<sub>2</sub>

**Disulfide Bridge:** Cys1 - Cys6

Standard  
Curve Image



Thursday 20 December, 2012

11983465 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.