

Secretin (Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-067-04	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Secretin

Normal Plasma Peptide/Protein Levels

Protocol	Assay Protocol (PDF)	
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction	
Catalog #	EK-067-04	
Sensitivity	Minimum Detectable Concentration = 0.22 ng/ml	
Range	0-100ng/ml	
Linear Range	0.22-2.36 ng/ml	
Cross Reactivity	Peptide	%Cross-reactivity
	Secretin (Mouse)	100
	Secretin (Rat)	100
	Secretin (Human)	100
	Secretin (Porcine)	100
	Secretin (Chicken)	0
	VIP (Human, Rat, Porcine)	0
	CCK (26-33) (Non Sulfated)	0
	Somatostatin	0
	NPY (Porcine)	0
	Met-Enkephalin	0
	Insulin (Human)	0
	Substance P	0
	PHI (Porcine)	0
	Amylin Amide (Human)	0
Sequence	His - Ser - Asp - Gly - Met - Phe - Thr - Ser - Glu - Leu - Ser - Arg - Leu - Gln - Asp - Ser - Ala - Arg - Leu - Gln - Arg - Leu - Leu - Gln - Gly - Leu - Val - NH ₂	

Standard
Curve Image

