

Secretin (Porcine, Bovine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-067-05	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

Secretin

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-067-05

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.18ng/ml

Range 0-100ng/ml

Linear Range 0.18-1.56ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Secretin (Porcine)	100
	Secretin (Human)	100
	Secretin (Rat)	100
	Secretin (Chicken)	0
	VIP (Human, Rat, Porcine)	0
	CCK (26-33) (Non Sulfated)	0
	Somatostatin	0
	NPY (Porcine)	0
	Met-Enkephalin	0
	Insulin (Human)	0
	Substance P	0
	PHI (Porcine)	0
	Amylin Amide (Human)	0

Sequence His - Ser - Asp - Gly - Thr - Phe - Thr - Ser - Glu - Leu - Ser - Arg - Leu - Arg - Asp - Ser - Ala - Arg - Leu - Gln - Arg - Leu - Leu - Gln - Gly - Leu - Val - NH₂

Standard
Curve Image

