

Irisin Precursor, C-terminal 48mer/ FNDC5 (165-212) (Human) / FNDC5 (162-209) (Rat, Mouse) - EIA kit

Catalog#	Standard size	Price	Order
EK-067-19	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

[Irisin - Antibody Stainings, Sequences & References](#)

[Irisin Synthesis and Its Applications](#)

Protocol	Assay Protocol (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	EK-067-19
Standard Size	96 Wells
Sensitivity	Minimum detection level: 0.1 ng/ml
Range	0-100ng/ml
Linear Range	0.21-2.01 ng/ml

Standard Curve Image

