

# CCK-33 (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-069-02	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

[All About Cholecystokinin](#)

[Normal Plasma Peptide/Protein Levels](#)

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-069-02

**Standard Size** 96 well

**Sensitivity** Minimum Detectable Concentration = 0.11 ng/ml

**Range** 0-100 ng/ml

**Linear Range** 0.11-2.33 ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	CCK-33 (Human)	100
	CCK-33 (Human) (Non-Sulfated)	42.2
	CCK (1-21) (Human)	42.2
	CCK-8 ( Human, Rat, Mouse) (26-33)	5.5
	CCK (26-33) (Non-Sulfated)	2.36
	CCK-8 (Non-Sulfated)	1
	CCK (1-21) (Porcine)	0
	Caerulein	0
	Gastrin-1 (Human)	0
	Pancreatic Polypeptide (Human)	0
	VIP (Human, Porcine, Rat)	0

**Sequence** Lys - Ala - Pro - Ser - Gly - Arg - Met - Ser - Ile - Val - Lys - Asn - Leu - Gln - Asn - Leu - Asp - Pro - Ser - His - Arg - Ile - Ser - Asp - Arg - Asp - Tyr(SO<sub>3</sub>H) - Met - Gly - Trp - Met - Asp - Phe - NH<sub>2</sub>

Standard  
Curve Image

