

CCK (26-33) (Non-Sulfated) (Human, Rat, Mouse, Canine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-069-04	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

[All About Cholecystokinin](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-069-04

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.06 ng/ml

Range 0-100ng/ml

Linear Range 0.06-0.85 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CCK (26-33) (Non-Sulfated) (Human, Rat, Mouse)	100
	CCK (Human) (1-21)	75
	CCK-33 (Porcine)	100
	CCK-33 (Sulfated) (Human)	42.9
	CCK-33 (Non-Sulfated) (Human)	77.8
	Caerulein	100
	Gastrin-1 (Human)	100
	Big Gastrin-1 (Human)	100
	CCK (26-33) (Sulfated)	3.9
	CCK (27-33)	100
	CCK (30-33)	12.8
	Pancreatic Polypeptide (Human)	0
	VIP (Human, Porcine, Rat)	0

Sequence Asp - Tyr - Met - Gly - Trp - Met - Asp - Phe - NH₂

