

Urotensin II (Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-071-08	1 kit	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

References

Urotensin II

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-071-08

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.11 ng/ml

Range 0-100 ng/ml

Linear Range 0.11-1.09 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Urotensin II (Human)	100
	Urotensin II (Rat)	100
	Urotensin II (Mouse)	100
	Endothelin-1	0
	Angiotensin II	0
	ADM (Rat)	0
	PAMP-20 (Human)	0
	α -ANP-28 (Human)	0
	CGRP (Rat)	0
	CNP-22	0
	BNP-25 (Rat)	0
	BNP-32 (Rat)	0
	Bradykinin	0

Sequence pGlu - His - Gly - Ala - Ala - Pro - Glu - Cys - Phe - Trp - Lys - Tyr - Cys - Ile

Disulfide Bridge: Cys8 - Cys13

