

TNF-Alpha, Recombinant (Human) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-072-28	1 kit	\$420.00	ORDER

Product description

More information

References

TNF (TUMOR NECROSIS FACTOR)

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\) *](#)

* - For Reference Only. Dilution conditions within this protocol may differ from each lot number. Please Contact us for more information on the protocol.

Sample Preparation [Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-072-28

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 43pg/ml

Linear Range 43-500 pg/ml

Cross Reactivity	peptide	%Cross-reactivity
	TNF- α (Human)	100
	TNF- α (Rat)	9.7
	TNF- α (Mouse)	7.2
	IL -6 (Human)	0
	Adiponectin (Human)	0
	Adiponectin (Mouse)	0
	Resistin (Human)	0
	Resistin (Mouse)	0
	Leptin (Human)	0
	Insulin (Human)	0
	C-Peptide (Human)	0
	GLP1 (Human)	0
	Amylin (Human)	0

