

Angiotensin I (Human, Rat, Mouse), Extraction Free EIA Kit

Catalog#	Standard size	Price	Order
EKE-002-01	1 kit	\$480.00	<input type="button" value="ORDER"/>

Product description

More information

Angiotensin I, II, III: Research Abstracts, Antibody Stainings, Kit Standard Curves, Sequences, etc.

Normal Plasma Peptide/Protein Levels

Sample Preparation	Blood collection Tissue extraction CSF extraction
Catalog #	EKE-002-01
Standard Size	96 Wells
Sensitivity	Minimum Detectable Concentration = 0.1ng/ml
Range	0-100 ng/ml
Linear Range	0.07-1.05 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Angiotensin I (Human, Rat, Mouse, Porcine, Canine, Guinea Pig)	100
	[Asn1, Val5, Asn9]-Angiotensin I (Salmon)	100
	[Val5, Asn9]-Angiotensin I (Bullfrog)	100
	[Val5, Ser9]-Angiotensin I (Fowl)	100
	Angiotensin II (Human, Rat, Mouse, Porcine, Canine)	0
	Angiotensin III (Human)	0
	Renin Substrate	0
	α-ANP (Human, Canine)	0
	BNP-32 (Human)	0
	[Arg8]-Vasopressin	0
	Endothelin-I (Human, Porcine, Rat, Mouse, Canine, Bovine)	0

Sequence Asp - Arg - Val - Tyr - Ile - His - Pro - Phe - His - Leu

Standard Curve Image



