

Angiotensin II (Human, Rat, Mouse) - Extraction Free EIA Kit

Catalog#	Standard size	Price	Order
EKE-002-12	1 kit	\$480.00	<input type="button" value="ORDER"/>

Product description

More information

Angiotensin I, II, III: Research Abstracts, Antibody Stainings, Kit Standard Curves, Sequences, etc.

Normal Plasma Peptide/Protein Levels

Sample Preparation
[Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EKE-002-12

Standard Size 96 Wells

Sensitivity Minimum Detectable Concentration = 0.08ng/ml

Range 0-100 ng/ml

Linear Range 0.08-1.55 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Angiotensin II (Human, Rat, Mouse, Porcine, Canine)	100
	[Val5]-Angiotensin II (Human)	100
	Angiotensin III (Human)	100
	Renin Substrate	0
	Angiotensin I (Human)	<1
	Angiotensin (1-7)	0
	Angiotensin I (1-9) (Human)	1
	Angiotensin II(1-7) (Human)	0
	α -ANP (1-28) (Human, Canine)	0
	[Arg8]-Vasopressin	0
	[Sar1, Ile8]-Angiotensin II	0

Sequence Asp - Arg - Val - Tyr - Ile - His - Pro - Phe

Standard Curve Image



