

# ANP, alpha (1-28) (Human), Extraction Free EIA Kit

Catalog#	Standard size	Price	Order
EKE-005-06	1 kit	\$480.00	<input type="button" value="ORDER"/>

## Product description

## More information

ANP (Atrial Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

**Sample Preparation**  
Blood collection  
Plasma extraction  
Tissue extraction  
CSF extraction

**Catalog #** EKE-005-06

**Standard Size** 96 Wells

**Sensitivity** Minimum Detectable Concentration = 0.10 ng/ml

**Range** 0-100ng/ml

**Linear Range** 0.10-1.51ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>% Cross-reactivity</b>
	$\alpha$ -ANP (1-28)	100
	$\beta$ -ANP (Human)	100
	ANP (Rat)	0
	BNP-26 (Porcine)	< 0.1
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	BNP-32 (Human)	0
	BNP-45 (Rat)	0
	CNP-22 (Human, Porcine, Rat)	0

**Sequence** Ser - Leu - Arg - Arg - Ser - Ser - Cys - Phe - Gly - Gly - Arg - Met - Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Leu - Gly - Cys - Asn - Ser - Phe - Arg - Tyr

## Image

