

BNP-32 (Human), Extraction Free EIA Kit

Catalog#	Standard size	Price	Order
EKE-011-03	1 kit	\$480.00	<input type="button" value="ORDER"/>

Product description

More information

Brain Natriuretic Peptides (BNP)

Normal Plasma Peptide/Protein Levels

Sample Preparation
[Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EKE-011-03

Standard Size 96 Wells

Sensitivity Minimum Detectable Concentration = 0.2ng/ml


Range 0-100 ng/ml

Linear Range 0.2-5.6 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	BNP-32 (Human)	100
	BNP-32 (Rat)	<0.04
	BNP-45 (Rat)	0
	BNP-45 (Mouse)	0
	BNP-26 (Porcine)	0
	CNP-22 (Human)	0
	α-ANP (1-28) (Human, Canine, Ovine, Porcine)	0
	Endothelin-1 (Human, Rat, Mouse, Porcine, Canine, Bovine)	0
	Angiotensin II (Human, Rat, Mouse, Porcine, Canine)	0
	VIP (Human, Rat, Porcine)	0

Sequence Asp - Arg - Val - Tyr - Ile - His - Pro - Phe

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.