

CCK (26-33)(non-Sulfated)(Human, Rat, Mouse), Extraction Free EIA Kit

Catalog#	Standard size	Price	Order
EKE-069-04	1 kit	\$480.00	<input type="button" value="ORDER"/>

[Product description](#) [More information](#)

[All About Cholecystokinin](#)

[Normal Plasma Peptide/Protein Levels](#)

Sample Preparation
[Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EKE-069-04

Standard Size 96 Wells

Sensitivity Minimum Detectable Concentration = 0.04ng/ml

Range 0-100ng/ml

Linear Range 0.04-0.45ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CCK (26-33) (Non-Sulfated) (Human, Rat, Mouse)	100
	CCK-33 (Porcine)	100
	Caerulein	100
	Gastrin-1 (Human)	100
	Big Gastrin-1 (Human)	100
	CCK (27-33)	100
	CCK-33 (Non-Sulfated) (Human)	77.8
	CCK (Human) (1-21)	75
	CCK-33 (Sulfated) (Human)	42.9
	CCK (30-33)	12.8
	CCK (26-33) (Sulfated)	3.9
	Pancreatic Polypeptide (Human)	0
	VIP (Human, Porcine, Rat)	0

Sequence Asp - Tyr - Met - Gly - Trp - Met - Asp - Phe - NH₂

Image



