

ACTH (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-001-01	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

ACTH

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-001-01

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 42.7 pg/ml

Range 0-10000 pg/ml

Linear Range 42.7-1828 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	ACTH (Human)	100
	ACTH (Rat)	100
	ACTH (1-24) (Human)	73
	ACTH (7-38) (Human)	31
	Beta-Endorphin (Human)	1
	CRF (Human, Rat)	0
	Alpha-MSH	0
	Met-Enkephalin	0
	Alpha-ANP (1-28) (Human, Canine)	0
	BNP-32 (Human)	0

Standard
Curve Image

