

ACTH (Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-001-21	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

ACTH

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-001-21

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 28.3pg/ml

Range 0-10000 pg/ml

Linear Range 28.3-400 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	ACTH (Rat, Mouse)	100
	ACTH (Human)	100
	ACTH (18-39) (Human)	100
	ACTH (1-24) (Human)	0.01
	ACTH (7-38) (Human)	0
	b-Endorphin (Rat)	0
	a-MSH	0
	LH-RH	0
	PACAP-38 (Ovine)	0

Sequence Ser - Tyr - Ser - Met - Glu - His - Phe - Arg - Trp - Gly - Lys - Pro - Val - Gly - Lys - Lys - Arg - Arg - Pro - Val - Lys - Val - Tyr - Pro - Asn - Val - Ala - Glu - Asn - Glu - Ser - Ala - Glu - Ala - Phe - Pro - Leu - Glu - Phe

