

Angiotensin I (Human, Rat, Mouse, Canine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-002-01	1 kit	\$460.00	<input type="button" value="ORDER"/>


- Product description**
- More information
- References

Angiotensin I, II, III: Research Abstracts, Antibody Stainings, Kit Standard Curves, Sequences, etc.
 Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-002-01
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 16 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Angiotensin I (Human, Mouse, Canine, Rat)	100
	[Asn ¹ , Val ⁵ , Asn ⁹]-Angiotensin I (Salmon)	100
	[Val ⁵ , Asn ⁹]-Angiotensin I (Bullfrog)	100
	[Val ⁵ , Ser ⁹]-Angiotensin I (Fowl)	100
	Angiotensin II (Human, Mouse, Canine, Rat)	0
	Angiotensin III (Human)	0
	Renin Substrate	0
	Alpha-ANP (Human, Canine)	0
	BNP-32 (Human)	0
	[Arg ⁸]-Vasopressin	0
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine)	0

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.