

Angiotensin II (Human, Rat, Mouse, Canine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-002-12	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Angiotensin I, II, III: Research Abstracts, Antibody Stainings, Kit Standard Curves, Sequences, etc.](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-002-12

Standard Size 96 well plate


Sensitivity Minimum Detectable Concentration = 13 pg/ml

Range 0-10000 pg/ml

Linear Range 13-240 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Angiotensin II (Human, Canine, Rat)	100
	[Val ⁵]-Angiotensin II (Human)	100
	Angiotensin III (Human)	100
	Renin Substrate	0
	Angiotensin I (Human)	<1
	Angiotensin (1-7)	0
	Angiotensin I (1-9) (Human)	1
	Angiotensin II (1-7) (Human)	0
	Alpha-ANP (1-28) (Human, Canine)	0
	[Arg ⁸]-Vasopressin	0
	[Sar ¹ , Ile ⁸]-Angiotensin II	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$