

Orexin A (Human, Rat, Mouse, Porcine, Ovine, Bovine, Monkey) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-003-30	1 kit	\$460.00	<input type="button" value="ORDER"/>


- Product description**
- More information
- References

Orexins and Obesity

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-003-30
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 42.5 pg/ml
Range	0-10000 pg/ml
Linear Range	42.5-480 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Orexin A (Human, Bovine, Mouse, Rat)	100
	Orexin A (16-33)	0
	Orexin B (Human)	0
	Agouti-Related Protein (83-132)-Amide (Human)	0
	Neuropeptide (Human, Rat)	0
	Alpha-MSH	0
	Leptin (Human)	0

Standard Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.