

AGRP (83-132)-Amide (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-003-53	1 kit	\$460.00	ORDER

Product description

More information

References

AGRP (Human)

Canine Beta-Defensin CBD103

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-003-53

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 11.2 pg/ml


Range 0-10000 pg/ml

Linear Range 11.2 □ 108 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Agouti-Related Protein (83-132)-NH ₂ (Human)	100
	Agouti-Related Protein Form C-NH ₂ (Human)	100
	Agouti-Related Protein (82-131)- NH ₂ (Mouse)	65
	Leptin (Human)	0
	Orexin A (Human, Mouse, Rat)	0
	Orexin B (Human)	0
	Neuropeptide Y (Human, Rat)	0
	α-MSH	0
	MCH (Human, Mouse, Rat)	0
	CGRP (Rat)	0
	CGRP (Human)	0

Sequence Ser - Ser - Arg - Arg - Cys - Val - Arg - Leu - His - Glu - Ser - Cys - Leu - Gly - Gln - Gln - Val - Pro - Cys - Cys - Asp - Pro - Cys - Ala - Thr - Cys - Tyr - Cys - Arg - Phe - Phe - Asn - Ala - Phe - Cys - Tyr - Cys - Arg - Lys - Leu - Gly - Thr - Ala - Met - Asn - Pro - Cys - Ser - Arg - Thr - NH₂

Standard
Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type Human Extracted* Plasma (Pooled)

Linear Range ** 13.9-141 pg/ml

ED50 28.05 pg/ml

Sample Level 2.25 pg/ml