

# ANP, alpha (1-28) (Human, Porcine, Ovine, Canine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-005-06	1 kit	\$460.00	<input type="button" value="ORDER"/>

- Product description**
- More information
- References

ANP (Atrial Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-005-06


**Standard Size** 96 well plate

**Sensitivity** Minimum Detectable Concentration = 14.8 pg/ml

**Linear Range** 14.8-338 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Alpha-ANP (1-28) (Human, Ovine, Canine)	100
	Beta-ANP (Human)	100
	ANP (Rat)	0
	BNP-26 (Porcine)	<0.1
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	BNP-32 (Human)	0
	BNP-45 (Rat)	0
	CNP-22 (Human, Porcine, Rat)	0

Standard Curve Image



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.