

# ANP (Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-005-24	1 kit	\$460.00	<input type="button" value="ORDER"/>

## Product description

## More information

ANP (Atrial Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-005-24

**Standard Size** 96 well


**Sensitivity** Minimum Detectable Concentration = 10.5 pg/ml

**Range** 0-10000 pg/ml

**Linear Range** 10.5-776 pg/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	ANP (1-28) (Rat, Mouse)	100
	$\alpha$ -ANP (8-33) (Rat)	100
	$\beta$ -ANP (1-28) (Human)	100
	Urodilatin	100
	$\alpha$ -ANP (1-28) (Human, Canine, Porcine, Bovine)	43
	BNP-26 (Porcine)	< 0.1
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	BNP-32 (Human)	0
	BNP-45 (Rat)	0
	CNP-22 (Human, Porcine, Rat)	0

Standard  
Curve Image



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$