

Bradykinin (Human, Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-009-01	1 kit	\$460.00	<input type="button" value="ORDER"/>

- Product description
- More information

- Bradykinin and Related Peptides
- Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)	
Sample Preparation	Blood collection	
	Plasma extraction	
	Tissue extraction	
	CSF extraction	
Catalog #	FEK-009-01	
Standard Size	96 well plate	
Sensitivity	Minimum Detectable Concentration = 47.5 pg/ml	
Range	0-10000 pg/ml	
Linear Range	47.5-1143 pg/ml	
Cross Reactivity	Peptide	%Cross-reactivity
	Bradykinin	100
	[Lys]-Bradykinin	100
	Biotinyl-Bradikinin	100
	[Tyr ^o]-Bradykinin	100
	[Des-Arg ⁷]-Bradykinin	0
	Bradykinin Potentiator B	0
	Angiotensin II (Human)	0
	Met-Enkephalin	0
	Substance P	0
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine)	0
	TRH	0

