

Adrenomedullin (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-010-01	1 kit	\$485.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Adrenomedullin](#)

[Normal Plasma Peptide/Protein Levels](#)

[Adrenomedullin-5 Porcine](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-010-01

Standard Size 96 well


Sensitivity Minimum Detectable Concentration = 31.6 pg/ml

Range 0-10000 pg/ml

Linear Range 31.6-428 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Adrenomedullin (Human)	100
	Adrenomedullin (Rat)	100
	Adrenomedullin (13-52) (Human)	100
	Adrenomedullin (1-50) (Mouse)	100
	Adrenomedullin-Gly (Human)	<1
	Amylin (Human)	0
	CGRP (Human)	0
	CGRP II (Human)	0
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine)	0
	a-ANP (1-28) (Human, Canine)	0
	BNP-32 (Human)	0
	CNP-22 (Human, Porcine, Rat)	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$