

BNP-32 (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-011-03	1 kit	\$460.00	<input type="button" value="ORDER"/>
Product description		More information	References

Brain Natriuretic Peptides (BNP)

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)


Catalog # FEK-011-03

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 31.8 pg/ml

Cross Reactivity	Peptide	% Crossreactivity
	BNP-32 (Human)	100
	BNP-32 (Rat)	<0.04
	BNP-45 (Rat)	0
	BNP-26 (Porcine)	0
	a-ANP (1-28) (Human, Canine)	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	Angiotensin II (Human)	0
	VIP (Human, Rat, Porcine)	0

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.