

# BNP-32 (Rat) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-011-14	1 kit	\$460.00	<input type="button" value="ORDER"/>

## Product description

## More information

Brain Natriuretic Peptides (BNP)

Normal Plasma Peptide/Protein Levels

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-011-14

**Standard Size** 96 well


**Sensitivity** Minimum Detectable Concentration = 61 pg/ml

**Range** 0-10000 pg/ml

**Linear Range** 61-1226 pg/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	BNP-32 (Rat)	100
	BNP-45 (Rat)	100
	BNP-32 (Human)	0
	BNP-26 (Porcine)	0
	CNP-53 (Porcine, Rat)	0
	ANP (Rat, Mouse)	0
	ET-1 (Human, Rat, Porcine, Canine)	0
	VIP (Human, Rat, Porcine)	0
	Angiotensin II (Human)	0

Standard  
Curve Image



Log-Logit Fit:  $y=(A-D)/(1+(x/C)^B)+D$