

BNP-45 (Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-011-23	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

Brain Natriuretic Peptides (BNP)

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-011-23

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 50 pg/ml


Range 0-10000 pg/ml

Linear Range 50-718 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	BNP-45 (Mouse)	100
	BNP-45 (Rat)	0
	BNP-32 (Porcine)	0
	BNP-26 (Porcine)	0
	BNP-32 (Human)	0
	BNP-32 (Canine)	0
	ANP (Rat)	0
	A-ANP (1-28) (Human, Canine)	0
	Endothelin-1	0
	Angiotensin II (Human)	0

Sequence Ser - Gln - Gly - Ser - Thr - Leu - Arg - Val - Gln - Gln - Arg - Pro - Gln - Asn - Ser - Lys - Val - Thr - His - Ile - Ser - Ser - Cys - Phe - Gly - His - Lys - Ile - Asp - Arg - Ile - Gly - Ser - Val - Ser - Arg - Leu - Gly - Cys - Asn - Ala - Leu - Lys - Leu - L

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$