

# Calcitonin (Salmon) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-014-09	1 kit	\$460.00	<input type="button" value="ORDER"/>

## Product description

## More information

Calcitonin Receptor-Stimulating Peptide (CRSP)

Normal Plasma Peptide/Protein Levels

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-014-09

**Standard Size** 96 well

**Sensitivity** Minimum Detectable Concentration = 59pg/ml

**Range** 0-10000 pg/ml

**Linear Range** 59-564pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Calcitonin (Salmon)	100
	Calcitonin (Eel)	100
	[Asu1,7]-Calcitonin (Eel)	100
	Calcitonin (Chicken)	100
	Calcitonin (Human)	0
	Calcitonin (Rat)	0
	Katacalcin (PDN-21)	0
	Calcitonin C-Terminal Adjacent Peptide (Rat)	0

**Sequence** Cys - Ser - Asn - Leu - Ser - Thr - Cys - Val - Leu - Gly - Lys - Leu - Ser - Gln - Glu - Leu - His - Lys - Leu - Gln - Thr - Tyr - Pro - Arg - Thr - Asn - Thr - Gly - Ser - Gly - Thr - Pro - NH<sub>2</sub> [Disulfide Bridge: Cys1 - Cys7]

Standard  
Curve Image

