

CGRP (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-015-02	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References


Calcitonin Gene Related Peptide (CGRP)

Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-015-02
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 60 pg/ml
Range	0-10000 pg/ml
Linear Range	60-1186 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CGRP (Human)	100
	CGRP (8-37) (Human)	100
	CGRP (Chicken)	100
	CGRP (Rat)	100
	CGRP II (Human)	100
	CGRP II (Rat)	100
	Calcitonin (Human)	0
	CAP (Rat)	0
	Amylin (Human)	0
	Amylin-Amide (Human)	0
	Amylin (Rat)	0
	Substance P	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$