


CGRP II (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-015-07	1 kit	\$460.00	<input type="button" value="ORDER"/>
Product description	More information	References	

Calcitonin Gene Related Peptide (CGRP)
Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)		
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction		
Catalog #	FEK-015-07		
Standard Size	96 well		
Sensitivity	Minimum Detectable Concentration = 29.2 pg/ml		
Range	0-10000 pg/ml		
Linear Range	29.2 668 pg/ml		
Cross Reactivity	Peptide	%Cross-reactivity	
	CGRP II (Human)	100	
	CGRP (8-37) (Human)	100	
	CGRP (Chicken)	100	
	CGRP (Rat)	100	
	CGRP (Human)	100	
	CGRP II (Rat)	100	
	Calcitonin (Human)	0	
	CAP (Rat)	0	
	Amylin (Human)	0	
	Amylin Amide (Human)	0	
	Amylin (Rat)	0	
	Substance P	0	
Sequence	Ala - Cys - Asn - Thr - Ala - Thr - Cys - Val - Thr - His - Arg - Leu - Ala - Gly - Leu - Leu - Ser - Arg - Ser - Gly - Gly - Met - Val - Lys - Ser - Asn - Phe - Val - Pro - Thr - Asn - Val - Gly - Ser - Lys - Ala - Phe - NH2 [Disulfide Bridge: Cys2 - Cys7]		

Standard Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type	Human Extracted* Plasma (Pooled)
Linear Range **	35.1-1140 pg/ml
ED50	125.62 pg/ml
<u>Sample Level</u>	31.5 pg/ml

Thursday 20 December, 2012

11984094 requests since Wednesday 05 April, 2006

Copyright © 2012 [PHOENIX PHARMACEUTICALS, INC.](#)