

CGRP (Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-015-09	1 kit	\$460.00	<input type="button" value="ORDER"/>


Product description **More information**

Calcitonin Gene Related Peptide (CGRP)
 Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-015-09
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 46 pg/ml
Range	0-10000 pg/ml
Linear Range	46-452 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CGRP (Rat)	100
	CGRP II (Rat)	78.6
	CGRP II (Human)	20.1
	CGRP (Human)	15.1
	Amylin-Amide (Rat)	<0.01
	Amylin (Human)	<0.001
	Calcitonin (Rat)	0
	Somatostatin-14	0

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.