

Corticotropin Releasing Factor (CRF) (Human, Rat, Mouse, Canine, Feline) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-019-06	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information


References

Corticotropin Releasing Factor (CRF)

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-019-06
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 45 pg/ml
Range	0-10000 pg/ml
Linear Range	45-868 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CRF (Human, Rat)	100
	CRF (Ovine)	0
	CRF (Bovine)	0
	Prepro-CRF (125-151) (Human)	0
	Sauvagine (Frog)	0
	PACAP-38 (Human, Rat, Ovine)	0
	LH-RH	0
	ACTH (Human)	0
	[Arg ⁸]-Vasopressin	0
	BNP-45 (Rat)	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$