

Endothelin-1 (ET-1) (Human, Rat, Mouse, Porcine, Bovine, Canine, Rabbit, Monkey) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-023-01	1 kit	\$460.00	ORDER

Product description

More information

References

[Endothelin and Endothelin-like Peptides](#)

[Preproendothelin-1: Placenta, Lung, Umbilical Vein \(Human\)](#)

[Endothelin-1 Precursor \(ET-1\) Sequence](#)


[Preproendothelin-1 \(Porcine\)](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-023-01
Standard Size	96 well
Sensitivity	Minimum Detectable Concentration = 12.9 pg/ml
Linear Range	17.8-185 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine, Rabbit, Monkey)	100
	Endothelin-2 (Human, Porcine, Rat, Mouse, Canine, Bovine)	100
	Endothelin-3 (Human, Porcine, Rat, Mouse, Canine, Bovine)	100
	Big Endothelin-1 (Human)	16
	Big Endothelin-2 (Human)	7
	Angiotensin II (Human)	0
	[Arg ⁸]-Vasopressin	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$