

Gastrin I (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-027-04	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

[Gastrin Precursor Schematic \(Human\)](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-027-04

Standard Size 96 well


Sensitivity Minimum Detectable Concentration = 7.5 pg/ml

Range 0-10000 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Gastrin I (Human)	100
	Big Gastrin I (Human)	100
	CCK (26-33) (Sulfated)	3
	CCK (26-33) (Non-Sulfated)	8
	Gastrin I (Chicken)	0.5
	Pancreatic Polypeptide (Human)	0
	GRP (Human)	0
	GIP (Porcine)	0
	Amylin Amide (Human)	0

Sequence pGlu - Gly - Pro - Trp - Leu - Glu - Glu - Glu - Glu - Glu - Ala - Tyr - Gly - Trp - Met - Asp - Phe - NH₂

Standard
Curve Image



Log-Logit Fit: $y=(A-D)/(1+(x/C)^B)+D$