

GLP-1 (7-36)-Amide (Human, Rat, Mouse, Porcine, Bovine, Canine, Ovine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-028-11	1 kit	\$395.00	<input type="button" value="ORDER"/>

Product description

More information

References

GLP-1 & GLP-2

Normal Plasma Peptide/Protein Levels

Neuron Function Modulation

Protocol [Fluorescent EIA Kit \(PDF\)](#)

*Please Note: For best results, use DPP-IV instead of Aprotinin.

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-028-11

Standard Size 96 well plate

Sensitivity Minimum Detectable Concentration = 27.6 pg/ml

Range 0-10000 pg/ml

Linear Range 27.6-676 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	GLP-1 (7-36) (Human, Rat, Mouse)	100
	GLP-1 (9-36) Amide	100
	GLP-1 (1-36) Amide (Human)	100
	GLP-1 (7-37) (Human, Rat, Mouse)	0.4
	GLP-2 (Human)	0
	Glucagon (Human, Rat, Mouse, Porcine, Bovine)	0
	CGRP	0
	GIP	0
	VIP	0
	Insulin	0
	Secretin	0
	Exendin-4	0

Standard
Curve Image

