

# GLP-2 (1-34) (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-028-14	1 kit	\$395.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

[GLP-1 & GLP-2](#)

[Normal Plasma Peptide/Protein Levels](#)

[Neuron Function Modulation](#)

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-028-14

**Standard Size** 96 wells

**Sensitivity** Minimum Detectable Concentration = 20 pg/ml


**Range** 0-10000 pg/ml

**Linear Range** 20-341 pg/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-Reactivity</b>
	GLP-2 (146-179) (Human)	100
	GLP-1 (7-36) (Human)	6
	GLP-1/ Pro-Glucagon (92-128) (Human, Rat)	5
	GLP-1 (7-37) (Human)	0
	Glucagon	0
	GIP	0
	VIP	0
	Insulin	0

**Sequence** His - Ala - Asp - Gly - Ser - Phe - Ser - Asp - Glu - Met - Asn - Thr - Ile - Leu - Asp - Asn - Leu - Ala - Ala - Arg - Asp - Phe - Ile - Asn - Trp - Leu - Ile - Gln - Thr - Lys - Ile - Thr - Asp - Arg

Standard  
Curve Image



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$