

Growth Hormone Releasing Factor (GHRF) (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-031-02	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Ghrelin Receptor Antagonist \[D-lys3\] / GHRP-6](#)

[Growth Hormone-Releasing Factor \(GHRF\) Prohormone Schematic](#)

[Normal Plasma Peptide/Protein Levels](#)

[JI-38 - Potent GHRH agonist](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)


Catalog # FEK-031-02

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 22.6 pg/ml

Linear Range 22.6-480 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Growth Hormone Releasing Factor (Human)	100
	Growth Hormone Releasing Factor (Porcine)	100
	Growth Hormone Releasing Factor (Bovine)	100
	Growth Hormone Releasing Factor (Ovine)	100
	Growth Hormone Releasing Factor (Rat)	0
	Growth Hormone Releasing Factor (1-40) (Human)	0
	Growth Hormone Releasing Factor (1-37) (Human)	0
	VIP (Human, Porcine, Rat)	0
	Growth Hormone (Human)	0
	Secretin (Human)	0
	GIP (Human)	0



Log-Logit Fit: $y=(A-D)/(1+(x/C)^B)+D$