

INSL 3 / RLF (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-035-27	1 kit	\$550.00	<button>ORDER</button>

- Product description
- More information

- Relaxin & Relaxin-Related Peptides
- Insulin-like 3 Peptide (INSL 3)
- Relaxin Evolution
- GSS: Relaxin Like Peptide

Protocol	Fluorescent EIA Kit (PDF)	
Sample Preparation	Blood collection	
	Plasma extraction	
	Tissue extraction	
	CSF extraction	
Catalog #	FEK-035-27	
Standard Size	96 well	
Sensitivity	Minimum Detectable Concentration = 19.6 pg/ml	
Range	0-10000 pg/ml	
Cross Reactivity	Peptide	%Cross-reactivity
	Insulin-Like 3 (INSL-3) (Human)	100
	Insulin-Like 3 (INSL-3) (Mouse)	0
	Insulin (Human)	0
	Insulin-Like 4 (INSL-4) (Human)	0
	Insulin-Like 5 (INSL-5) (Human)	0
	Insulin-Like 6 (INSL-6) (Human)	0
	Insulin-Like 7 (INSL-7) (Human)	0
	C-peptide (Human)	0
	Amylin (Human)	0
	Glucagon (Human)	0
	GLP-1 (7-37) (Human)	0
	GLP-2 (Human)	0
	Ghrelin (Human)	0
	LH-RH (Human, Porcine, Rat)	0
	CGRP (Human)	0
	AGRP (83-132) –Amide (Human)	0

