

MSH, Alpha (Human, Rat, Mouse) - Fluorescent EIA Kit


Catalog#	Standard size	Price	Order
FEK-043-01	1 kit	\$460.00	<input type="button" value="ORDER"/>
Product description		More information	References

Alpha-MSH (Alpha-Melanocyte Stimulating Hormone)

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-043-01
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 20 pg/ml
Range	0-10000 pg/ml
Linear Range	20-384 pg/ml

Cross Reactivity	Peptide	Crossreactivity
	alpha-MSH (Human, Rat, Mouse)	100
	Des Ac-alpha-MSH	79
	alpha-MSH (4-13) (Ac-(Cys4, d-Phe7, Cys10))	16.3
	Neuropeptide Y (Human, Rat)	0
	CART (55-102) (Human)	0
	AGRP (83-132)-NH2 (Human)	0
	ACTH (Human)	0
	Leptin (Human)	0

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.