

# Neurokinin A / Neuromedin L / Substance K (Human, Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-046-15	1 kit	\$460.00	<input type="button" value="ORDER"/>


## Product description

## References

<b>Protocol</b>	Fluorescent EIA Kit (PDF)
<b>Sample Preparation</b>	<a href="#">Blood collection</a> <a href="#">Plasma extraction</a> <a href="#">Tissue extraction</a> <a href="#">CSF extraction</a>
<b>Catalog #</b>	FEK-046-15
<b>Standard Size</b>	96 well plate
<b>Sensitivity</b>	Minimum Detectable Concentration = 26.6 pg/ml
<b>Range</b>	0-10000 pg/ml
<b>Linear Range</b>	26.6-577 pg/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	Neurokinin A	100
	Kassinin K	100
	Neuropeptide (Porcine)	100
	Neurokinin B	57
	Substance P	35
	Physalaemin	0.02
	Neuromedin B	0
	Neuropeptide Y (Porcine)	0
	ACTH (Human)	0
	[Arg <sup>8</sup> ]-Vaspressin	0

Standard  
Curve Image



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$