

# Neurokinin B / Neuromedin K (Human, Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-046-26	1 kit	\$460.00	<input type="button" value="ORDER"/>
<b>Product description</b>		<b>More information</b>	<b>References</b>

Neurokinin B {NKB} (Neuromedin K) Precursor Schematic (Human)

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-046-26

**Standard Size** 96 wells


**Sensitivity** Minimum Detectable Concentration = 24.4 pg/ml

**Range** 0-10000 pg/ml

**Linear Range** 24.4-608 pg/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>% Cross-reactivity</b>
	Neurokinin B	100
	Physalaemin	3.84
	Eledoisin	0
	Kassinin	0
	Substance P	0
	Neuromedin B (Porcine)	0
	Neuropeptide K (Porcine)	0

**Standard Curve Image**



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.