

Neurotensin (Human, Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-048-03	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

Neurotensin

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-048-03

Standard Size 96 well


Sensitivity Minimum Detectable Concentration = 57 pg/ml

Range 0-10000 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Neurotensin	100
	Neurotensin (Guinea Pig)	26
	Neuromedin N (Porcine)	0.03
	Substance P	0
	Kinetensin	0
	Bombesin	0

Sequence pGlu - Leu - Tyr - Glu - Asn - Lys - Pro - Arg - Arg - Pro - Tyr - Ile - Leu

Standard
Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type Human Extracted* Plasma (Pooled)

Linear Range ** 102-1063 pg/ml

ED50 331.5 pg/ml

Sample Level **39 pg/ml**