

NPY (Human, Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-049-03	1 kit	\$460.00	ORDER

Product description

More information

References

[Normal Plasma Peptide/Protein Levels](#)

[Neuropeptide Y \(NPY\)](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-049-03

Standard Size 96 well plate

Sensitivity Minimum Detectable Concentration = 11.9 pg/ml

Range 0-10000 pg/ml

Linear Range 11.9-531 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Neuropeptide Y (Human, Rat)	100
	Neuropeptide Y (Porcine)	100
	Neuropeptide Y (3-36) (Human, Rat, Mouse)	14.3
	Peptide YY (Human)	0
	Peptide YY (Porcine)	0
	Pancreatic Polypeptide (Human)	0
	VIP (Human, Rat, Porcine)	0
	Insulin (Human)	0
	Amylin (Human)	0
	Somatostatin	0

Standard
Curve Image

