

Chromogranin A, Prepro-(342-355) (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-053-26	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

Chromogranin A

Serpinin

Protocol	Fluorescent EIA Kit (PDF)										
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction										
Catalog #	FEK-053-26										
Standard Size	96 well										
Sensitivity	Minimum Detectable Concentration = 17.4 pg/ml										
Range	0-10000 pg/ml										
Linear Range	17.4-290 pg/ml										
Cross Reactivity	<table><thead><tr><th><i>Peptide</i></th><th><i>%Cross-reactivity</i></th></tr></thead><tbody><tr><td>Prepro-Chromogranin A (342-355) (Human)</td><td>100</td></tr><tr><td>Prepro-Chromogranin A (393-426) (Human)</td><td>0</td></tr><tr><td>Prepro-Chromogranin A (429-457) (Human)</td><td>0</td></tr><tr><td>Prepro-Chromogranin A (392-402) (Mouse)</td><td>0</td></tr></tbody></table>	<i>Peptide</i>	<i>%Cross-reactivity</i>	Prepro-Chromogranin A (342-355) (Human)	100	Prepro-Chromogranin A (393-426) (Human)	0	Prepro-Chromogranin A (429-457) (Human)	0	Prepro-Chromogranin A (392-402) (Mouse)	0
<i>Peptide</i>	<i>%Cross-reactivity</i>										
Prepro-Chromogranin A (342-355) (Human)	100										
Prepro-Chromogranin A (393-426) (Human)	0										
Prepro-Chromogranin A (429-457) (Human)	0										
Prepro-Chromogranin A (392-402) (Mouse)	0										
Sequence	Trp - Ser - Lys - Met - Asp - Gln - Leu - Ala - Lys - Glu - Leu - Thr - Ala - Glu										

Thursday 20 December, 2012

11984136 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.