

Serpinin (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-053-45	1 kit	\$460.00	<input type="button" value="ORDER"/>

[Product description](#) [More information](#)

[Chromogranin A](#)
[Serpinin](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation
[Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-053-45

Standard Size 96 Well EIA Kit


Sensitivity Minimum Detectable Concentration = 30.8 pg/ml

Range 0.01-10000 pg/ml

Linear Range 30.8-366 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Serpinin / Prepro-Chromogranin A (429-454) (Human)	100
	Serpinin-RRG / Prepro-Chromogranin A (429-457) (Human)	100
	Serpinin / Prepro-Chromogranin A (435-460) (Rat, Mouse)	0
	[pGlu]-Serpinin (Human, Rat, Mouse)	0

Standard Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$

