

Parathyroid Hormone (PTH) (1-34) (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-055-08	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

Parathyroid Hormone and its Related Protein

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-055-08

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 18.3 pg/ml


Range 0-10000 pg/ml

Linear Range 18.3-493 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	PTH (1-34) (Human)	100
	PTH (1-34) (Rat)	85.7
	VIP	0
	CGRP (Human)	0
	Calcitonin (Salmon)	0
	GLP-2	0
	PTH (Human) (53-84)	0

Sequence Ser - Val - Ser - Glu - Ile - Gln - Leu - Met - His - Asn - Leu - Gly - Lys - His - Leu - Asn - Ser - Met - Glu - Arg - Val - Glu - Trp - Leu - Arg - Lys - Lys - Leu - Gln - Asp - Val - His - Asn - Phe

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$