

Parathyroid Hormone-related Protein (PTH-RP) (1-34) (Human, Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-056-04	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

Parathyroid Hormone and its Related Protein

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-056-04

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 57 pg/ml


Range 0-10000 pg/ml

Linear Range 57-1031 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	PTH-rP (1-34) (Human)	100
	PTH-rP (1-37) (Human)	100
	PTH-rP (7-34) Amide (Human)	100
	PTH-rP (Chicken)	0
	PTH (1-34) (Human)	0
	PTH (1-34) (Rat)	0
	[Asn ⁷⁶]-PTH 184 (Human)	0
	PTH-rP (38-64) Amide (Human)	0
	PTH-rP (67-86) Amide (Human)	0
	PTH-rP (107-138) (Human)	0
	PTH-rP (109-141) (Human)	0

Sequence Ala - Val - Ser - Glu - His - Gln - Leu - Leu - His - Asp - Lys - Gly - Lys - Ser - Ile - Gln - Asp - Leu - Arg - Arg - Arg - Phe - Phe - Leu - His - His - Leu - Ile - Ala - Glu - Ile - His - Thr - Ala

Standard
Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type Human Extracted* Plasma (Pooled)

Linear Range ** 57-1031 pg/ml

ED50 178.87 pg/ml

Sample Level **29.1 pg/ml**