

Apelin-12 (Human, Rat, Mouse, Bovine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-057-23	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References


[Apelin - Research Abstracts, Antibody Stainings, Kits Standard Curve, Sequences, etc.](#)

Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-057-23
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 15.8 pg/ml
Range	0-10000 pg/ml
Linear Range	15.8-419 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Apelin-12 (Human)	100
	Apelin-13 (Human)	100
	Apelin-13 (Rat)	100
	Apelin-36 (Human)	100
	Apelin-36 (Rat)	100
	Alpha-ANP (Human)	0
	ANP (25-56) (N-ANP)	0
	ADM-52 (Human)	0
	BNP-32 (Human)	0
	CNP-22	0
	Ghrelin (Human)	0
	ET-1 (Human)	0
	Bradykinin	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$